

# NEOTEMP 1965

## Neonatal Skin Temperature Sensor



### For use with Dräger Medical Caleo Incubator and Babytherm 8004 / 8010 Radiant Warmer

(Replaces Dräger Medical  
Thermotrace Core Temperature Probe "Yellow"  
Part# MX11000)

#### Temp Accuracy

NovaMed baby-sensitive thermistor sensor guarantees temperature reading sensitivity with a stable, highly accurate measurement.

Molded plug assembly provides direct connection to the Caleo Incubator / Babytherm Warmer.

#### Direct Connection

Flexible 62" lead wire facilitates sensor placement on the infant, eliminating the need for interface/interconnect cable and ensuring precision and safety.

#### Ease of Use

Packaged individually in a Tyvek pouch with an insulated adhesive Mylar temperature probe cover (one-inch circular disk).

#### Gentle

Neonatal sensitive hypoallergenic adhesive/reflective Mylar probe cover provides insulation from heat and external sources ensuring precise temperature monitoring capabilities.



### NEOTEMP Neonatal Skin Temperature Sensor Disposable with Insulated Mylar Disc

- Designed for Dräger Medical Caleo / Babytherm 8004/8010.
- Replaces Dräger Medical MX-11000 (Thermotrace Core Temperature Probe "Yellow") and MX11002 (ThermoPad Insulated Adhesive).
- NEOTEMP 1965; Packaged 20/box, Order Ref. #10-1965-020  
NEOTEMP 1965; Packaged 60/case, Order Ref. #10-1965-060