

**NOVA***lite*  
**Fiber Optic Laryngoscopes**

**XENON Technology – Green System**



**Optimal Illumination – Reusable and Single Use**

---

**NOVAMEDUSA**

800-425-3535  
[www.NovaMed-USA.com](http://www.NovaMed-USA.com)



## Fiber Optic Laryngoscopes Ultra-Brite XENON Technology



### Illumination



Engineered with XENON light technology to facilitate superior illumination. Designed with an enlarged 4mm fiber optic bundle for improved light output.

### Compatibility



Utilizes the universal green system to ensure interchangeability; compatible with all green system handles and blades from competitive brands.

### Patient Safety



Satin finish blades enhance the ease of intubations and improve patient safety. Ensures enhanced patient comfort.

### Flexibility



Designed with an adjustable and entirely removable Fiber Optic bundle for enhanced cleaning benefits for improved safety and enhanced hospital risk management.

### Extensive Selection



Available in a comprehensive range of blade styles and sizes including Macintosh 3.5, Miller 1.5, and Wis-Hipple 1.5.



Fiber Optic Green System Handles  
Shown above (left to right)  
Medium, Penlight and Stubby

Reusable Handles Fiber Optic Green System		
Size	Batteries	Order #
Medium/ Adult	2C	14 - 600FO - 01
Penlight/ Pediatric	2AA	14 - 600FO - 02
Stubby	2AA	14 - 600FO - 03



Xenon Light Technology Replacement Bulbs	
Replacement Xenon Bulbs	Order #
Ultra Brite Xenon Bulbs 2.7V (6 pack)	14 - 700FO - 02 - 06



## NOVALITE Fiber Optic Reusable Laryngoscope Blades

Stainless Steel, Satin Finish – Green System

High quality stainless steel ensures long lasting performance.

Featuring an enlarged fiber optic bundle for optimal illumination.

### MACINTOSH FIBER OPTIC REUSABLE

Order #	Size / Description
14 - 100FO - MAC - 0	0, Neonate
14 - 100FO - MAC - 1	1, Infant
14 - 100FO - MAC - 2	2, Child
14 - 100FO - MAC - 3	3, Medium Adult
14 - 100FO - MAC - 3.5	3.5, Medium Adult +
14 - 100FO - MAC - 4	4, Adult
14 - 100FO - MAC - 5	5, Large Adult

### MILLER FIBER OPTIC REUSABLE

Order #	Size / Description
14 - 300FO - MIL - 00	00, Preemie
14 - 300FO - MIL - 0	0, Neonate
14 - 300FO - MIL - 1	1, Infant
14 - 300FO - MIL - 1.5	1.5, Infant +
14 - 300FO - MIL - 2	2, Child
14 - 300FO - MIL - 3	3, Medium Adult
14 - 300FO - MIL - 4	4, Adult

### FLIP TIP MACINTOSH FIBER OPTIC REUSABLE

Designed to facilitate difficult intubations with flexible tip adjustable through 70°

Order #	Size / Description
14 - 200FOFT - MAC - 2	2, Child
14 - 200FOFT - MAC - 3	3, Medium Adult
14 - 200FOFT - MAC - 4	4, Adult

### SPECIALTY BLADES FIBER OPTIC REUSABLE

	Order #	Size
Robert Shaw	14 - 400FO - RS - 0	0, Neonate
	14 - 400FO - RS - 1	1, Infant
Wis - Hipple	14 - 800FO - WH - 1.5	1.5, Infant +
Cranwall	14 - 400FO - CW - 2	2, Child
	14 - 400FO - CW - 3	3, Medium Adult
	14 - 400FO - CW - 4	4, Adult



## Fiber Optic Single Use Hybrid Laryngoscope Blade Stainless Steel Green System with Acrylic FO Bundle



- For safer intubation eliminating cross contamination.
- Affordably priced to eliminate cleaning requirements.
- Single use stainless steel blade with polished Acrylic bundle engineered to maximize illumination.
- Ideal for crash carts / intubation boxes, NICU, PICU, L&D and ER.
- Compatible with Green System Fiber Optic Handles.



<b>MACINTOSH FIBER OPTIC SINGLE USE</b>	
<b>Order # (10 pcs / Package)</b>	<b>Size / Description</b>
16 - 100FOMD - MAC - 0 - 010	0, Neonate
16 - 100FOMD - MAC - 1 - 010	1, Infant
16 - 100FOMD - MAC - 2 - 010	2, Child
16 - 100FOMD - MAC - 3 - 010	3, Medium Adult
16 - 100FOMD - MAC - 4 - 010	4, Adult
16 - 100FOMD - MAC - 5 - 010	5, Large Adult
<b>OXFORD FIBER OPTIC SINGLE USE</b>	
10 - 400FOMD - OX - 010	1, Infant

<b>MILLER FIBER OPTIC SINGLE USE</b>	
<b>Order # (10 pcs / Package)</b>	<b>Size / Description</b>
16 - 300FOMD - MIL - 00 - 010	00, Premie
16 - 300FOMD - MIL - 0 - 010	0, Neonate
16 - 300FOMD - MIL - 1 - 010	1, Infant
16 - 300FOMD - MIL - 2 - 010	2, Child
16 - 300FOMD - MIL - 3 - 010	3, Medium Adult
16 - 300FOMD - MIL - 4 - 010	4, Adult