

NEOTEMP 1950

Neonatal Skin Temperature Sensor



For use with
**GE/ Ohmeda Medical
GIRAFFE & GIRAFFE Omni-Bed
Incubators and
PANDA Warmers**

(Replaces GE/ Ohmeda Medical GIRAFFE Omni-Bed
#6600-0873-700 (pkg of 10) &
#6600-0874-700(pkg of 50))

Temp Accuracy

NovaMed baby-sensitive thermistor sensor guarantees temperature reading sensitivity with a stable, highly accurate measurement.

Molded plug assembly provides direct connection to the GIRAFFE incubator/radiant warmer.

Sensor Safety

Flexible lead wire facilitates sensor placement on the infant eliminating the need for interface/interconnect cable, ensuring precision and safety.


Ease of Use

Packaged individually in a Tyvek pouch with an insulated adhesive Mylar temperature probe cover (one inch circular disk).

Gentle

Neonatal sensitive hypoallergenic adhesive/reflective Mylar probe cover provides insulation from heat and external sources while further ensuring precise temperature monitoring capabilities.

NEOTEMP Neonatal Skin Temperature Sensor Disposable with Insulated Mylar Disc



- Designed exclusively for Ohmeda Medical Omni-Bed Giraffe & Giraffe Incubators and Giraffe/ Panda Radiant Warmers.
- Replaces Ohmeda Medical Omni-Bed Giraffe Disposable Skin Temperature Probe # 6600-0873-700 & # 6600-0874-700 (Reusable # 6600-0875-700).
- NEOTEMP 1950; Packaged 20/box, Order Ref. #10-1950-020
NEOTEMP 1950; Packaged 60/case, Order Ref. #10-1950-060.
- 3 Pin Din molded connector to fit the GE/Ohmeda Omni-Bed Giraffe /Giraffe /Panda.