


NOVALITE

LED Standard Laryngoscopes

The NEW Standard for Improved Illumination

- The ONLY LED SS blades available in Reusable & Single Use
- Whiter, brighter light for optimal task illumination
- 5X light output of standard incandescent blades
- Cool to the touch and low profile, ideal for neonatal



LED technology ensures whiter, brighter illumination with fewer shadows. The light emitting diode generates an increased number of photons which bounce off the inner surface in a mirrored effect. The result is increased light output with optimal, focused task illumination.

- Reduced energy consumption prolongs battery life
 - Compatible with Standard Conventional handles
 - Upgrade to improved illumination at a cost savings
-

LED BRITE Reusable Standard SS Laryngoscope Blades



- Superior Illumination** Engineered with a whiter, brighter light.
5X light output for optimal, focused task illumination.
Ensures outstanding light transmission and improved visualization.
- Compatibility** Compatible with Standard Conventional Handles.
- Improved Safety** LED remains cool at all times.
Directional light is optimal for task performance.
Satin finish blades ensure ease of intubations to enhance patient safety.
- Cost Effective** No light filaments are utilized, product efficiency is optimized.
Reduced energy consumption minimizes battery replacement.
High quality stainless steel ensures long lasting performance.
- Comprehensive Selection** Available in both MacIntosh & Miller styles in a complete range of sizes.

MACINTOSH LED BRITE STANDARD BLADES			
Reusable Stainless Steel, Satin Finish			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 6010R	MAC 0, Neonate	55	9
18 - 6011R	MAC 1, Infant	70	11
18 - 6012R	MAC 2, Child	78	13
18 - 6013R	MAC 3, Medium Adult	108	15
18 - 6014R	MAC 4, Adult	133	15
18 - 6015R	MAC 5, Large Adult	150	15

MILLER LED BRITE STANDARD BLADES			
Reusable Stainless Steel, Satin Finish			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 60200R	MIL 00, Premature	45	9.5
18 - 6020R	MIL 0, Neonate	55	9.5
18 - 6021R	MIL 1, Infant	80	9.5
18 - 6022R	MIL 2, Child	130	13
18 - 6023R	MIL 3, Medium Adult	170	13
18 - 6024R	MIL 4, Adult	180	13

LED BRITE Single Use Standard SS Laryngoscope Blades



Single Use SS

Whiter Brighter Light Engineered with a whiter, brighter light.
Improved directional, focused light with 5X light output.
Superior intubation conditions than with incandescent light sources.

Compatibility Compatible with Standard Conventional Handles.

Safe Prevents cross-infection.
Ideal for ER's or crash carts.
Individually wrapped, ready for use.
LED bulb remains cool at all times.
Directional light improves ease of intubations and patient safety.

Cost Effective Eliminates reprocessing/cleaning costs.
Reduced energy consumption minimizes battery replacement.

Comprehensive Selection Available in both MacIntosh & Miller styles in a complete range of sizes.

MACINTOSH LED BRITE STANDARD BLADES

Single Use Stainless Steel, Satin Finish

Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 7010D	MAC 0, Neonate	55	9
18 - 7011D	MAC 1, Infant	70	11
18 - 7012D	MAC 2, Child	78	13
18 - 7013D	MAC 3, Medium Adult	108	15
18 - 7014D	MAC 4, Adult	133	15
18 - 7015D	MAC 5, Large Adult	150	15

MILLER LED BRITE STANDARD BLADES

Single Use Stainless Steel, Satin Finish

Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 7020D	MIL 00, Premature	45	9.5
18 - 7020D	MIL 0, Neonate	55	9.5
18 - 7021D	MIL 1, Infant	80	9.5
18 - 7022D	MIL 2, Child	130	13
18 - 7023D	MIL 3, Medium Adult	170	13
18 - 7024D	MIL 4, Adult	180	13

To Order Replacement or Additional NOVALITE Standard / Conventional Laryngoscopes Handles

Shown from left to right:
Medium / Adult Handle
Penlight / Pediatric Handle
Stubby Handle
Large Handle



- High quality stainless steel construction with polished finish.
- Compatible with other standard / conventional blades.
- Available in a complete range of sizes as listed below.

Order #	Description	Batteries
14 - 600S - 01	Medium Adult Standard Handle	2 C
14 - 600S - 02	Penlight/Pediatric Standard Handle	2 AA
14 - 600S - 03	Stubby Standard Handle	2 AA
14 - 600S - 06	Large Standard Handle	2 D

Always Innovative • Always Reliable • Always Trusted

NOVAMEDUSA

800-425-3535
www.NovaMed-USA.com