


# NOVA*lite*

## LED Standard Laryngoscopes

The NEW Standard for Improved Illumination

- The ONLY LED SS blades available in Reusable & Single Use
- Whiter, brighter light for optimal task illumination
- 5X light output of standard incandescent blades
- Cool to the touch and low profile, ideal for neonatal



LED technology ensures whiter, brighter illumination with fewer shadows. The light emitting diode generates an increased number of photons which bounce off the inner surface in a mirrored effect. The result is increased light output with optimal, focused task illumination.

- Reduced energy consumption prolongs battery life
- Compatible with Standard Conventional handles
- Upgrade to improved illumination at a cost savings



# LED BRITE Reusable Standard SS Laryngoscope Blades



- Superior Illumination** Engineered with a whiter, brighter light.  
 5X light output for optimal, focused task illumination.  
 Ensures outstanding light transmission and improved visualization.
- Compatibility** Compatible with Standard Conventional Handles.
- Improved Safety** LED remains cool at all times.  
 Directional light is optimal for task performance.  
 Satin finish blades ensure ease of intubations to enhance patient safety.
- Cost Effective** No light filaments are utilized, product efficiency is optimized.  
 Reduced energy consumption minimizes battery replacement.  
 High quality stainless steel ensures long lasting performance.
- Comprehensive Selection** Available in Macintosh 0-5 & Miller 00-4. Also available in specialty blades, including Wis Hipple 1.5 and Robert Shaw 0-1.

<b>MACINTOSH LED BRITE STANDARD BLADES</b>			
<b>Reusable Stainless Steel, Satin Finish</b>			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 6010R	MAC 0, Neonate	55	9
18 - 6011R	MAC 1, Infant	70	11
18 - 6012R	MAC 2, Child	88	13
18 - 6013R	MAC 3, Medium Adult	108	15
18-6013.5R	Mac 3.5, Med. Adult +	115	15
18 - 6014R	MAC 4, Adult	133	15
18 - 6015R	MAC 5, Large Adult	150	15

<b>MILLER LED BRITE STANDARD BLADES</b>			
<b>Reusable Stainless Steel, Satin Finish</b>			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 60200R	MIL 00, Premature	45	9.5
18 - 6020R	MIL 0, Neonate	55	9.5
18 - 6021R	MIL 1, Infant	80	9.5
18-6021.5R	MIL 1.5, Infant +	90	9.5
18 - 6022R	MIL 2, Child	130	13
18 - 6023R	MIL 3, Medium Adult	170	13
18 - 6024R	MIL 4, Adult	180	13



# LED BRITE Single Use Standard SS Laryngoscope Blades



Single Use  
100% Stainless Steel

- Whiter Brighter Light**      Engineered with a whiter, brighter light.  
Delivers 5X light output of standard incandescent blades.
- Compatibility**              Compatible with Standard Conventional Handles.
- 100% Stainless**              The ONLY 100% Stainless Steel Single Use LED available.
- Safe**                              Prevents cross-contamination – Improves reliability.  
Bulb is fixed to blade for improved patient safety.  
LED is cool in operation – safer for neonates.
- Cost Effective**                Eliminates reprocessing/cleaning costs.  
Reduced energy consumption minimizes battery replacement.
- Comprehensive Selection**    Available in both MacIntosh 0-5 & Miller 00-4.

<b>MACINTOSH LED BRITE STANDARD BLADES</b> Single Use Stainless Steel, Satin Finish			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 7010D	MAC 0, Neonate	55	9
18 - 7011D	MAC 1, Infant	70	11
18 - 7012D	MAC 2, Child	78	13
18 - 7013D	MAC 3, Med Adult	108	15
18 - 7013.5D	MAC 3.5, Med Adult +	115	15
18 - 7014D	MAC 4, Adult	133	15
18 - 7015D	MAC 5, Large Adult	150	15

<b>MILLER LED BRITE STANDARD BLADES</b> Single Use Stainless Steel, Satin Finish			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 70200D	MIL 00, Premature	45	9.5
18 - 7020D	MIL 0, Neonate	55	9.5
18 - 7021D	MIL 1, Infant	80	9.5
18 - 7021.5D	MIL 1.5, Infant +	90	9.5
18 - 7022D	MIL 2, Child	130	13
18 - 7023D	MIL 3, Medium Adult	170	13
18 - 7024D	MIL 4, Adult	180	13

To Order Replacement or Additional

**NOVA**lite  
Standard / Conventional  
Laryngoscopes Handles



Shown from left to right:  
Medium / Adult Handle  
Penlight / Pediatric Handle  
Stubby Handle  
Large Handle

- High quality stainless steel construction with polished finish.
- Compatible with other standard / conventional blades.
- Available in a complete range of sizes as listed below.

Order #	Description	Batteries
14 - 600S - 01	Medium Adult Standard Handle	2 C
14 - 600S - 02	Penlight/Pediatric Standard Handle	2 AA
14 - 600S - 03	Stubby Standard Handle	2 AA
14 - 600S - 06	Large Standard Handle	2 D

**Light Emitting Diode (LED) Replacement Lamps:**

Order #	Description	Packaging
14 - 700LED - 06	Standard Reusable LED Replacement Lamp	6 per case

Always Innovative • Always Reliable • Always Trusted