

NEOTEMP

Neonatal Skin Temperature Sensors
for Precise Temperature Measurement
With Baby Sensitive Accuracy




Designed to fit Incubators and Warmers from
GE/Ohmeda Medical,
Dräger Medical/Air Shields,
ATOM Medical
International Biomedical/AirBorne

NOVAMEDUSA

914-789-2100 / 800-425-3535
www.NovaMed-USA.com

NEOTEMP

Neonatal Skin Temperature Sensors



- Baby-sensitive thermistor sensor guarantees temperature reading sensitivity with a stable, highly accurate measurement.
- Comprehensive range of temperature sensors designed for GE / Ohmeda, Dräger / Air Shields, ATOM Medical and International Biomedical transport equipment.
- Molded plug assembly provides direct connection to the incubator / warmer, eliminating the need for an interface cable.
- Packaged with a KOALA disk cover for direct placement on infant, ensuring precise temperature reading.



KOALA 1" Round Adhesive Replacement Disk Covers

- Hypoallergenic adhesive cover.
- Reflective Mylar provides insulation from heat and external sources.

Product #	10 - 1900N - 300
Packaging	300 / Box

For GE / Ohmeda Medical Incubators & Radiant Warmers



NEOTEMP 1930 Skin Temperature Sensor

For use with Ohmeda Care Plus, IC, CG Incubators & Infant Warming Systems (3000, 3300, 3500 & 5000)

Replaces OEM #LA003, #LA005, #6600-0196-700 & #6600-0208-700



NEOTEMP 1950 Skin Temperature Sensor

For use with GE GIRAFFE and GIRAFFE Omni-Bed Incubators & PANDA Warmers

Replaces OEM #6600-0873-700 & #6600-0874-700

For Dräger Medical / AirShields Incubators & Radiant Warmers



NEOTEMP 1910 Skin Temperature Sensor

For use with ISOLETTE Incubators (C-100, C-450, C-550)
& PM-78 Radiant Warmers

Replaces OEM Premie Probe 1 #68-209-20/MU06933 and
#68-209-30/MU06935



NEOTEMP 1960 Skin Temperature Sensor

For use with ISOLETTE C2000 & 8000 Incubators &
VERSALET 770 Care Center

Replaces OEM Premie Probe 5/5-L #83-104-15/MU12551



NEOTEMP 1965 Skin Temperature Sensor

For use with CALEO Incubators & Babytherm Warmers 8004/8010

Replaces OEM #MX11000



NEOTEMP 1970 Skin Temperature Sensor

For use with RESUSCITAIRE Radiant Warmers

Replaces OEM Premie Probe 3 #81-300-006,
#81-300-008/MU11228 and #81-300-009/MU11229

For ATOM Medical Incubators & Radiant Warmers



NEOTEMP 2000D (Single Use) or 2000R (Reusable) Skin Temperature Sensors

For use with Atom Medical Incubators & Radiant Warmers

Replaces OEM NF-023, NF-024, #60882, #60883 &
#60084 with connector

NOVATEMP


**Skin Temperature Sensors & Temp Adapter Cable
For International Biomedical/AirBorne and
Dräger Transport Incubators**

Order #	Product	Packaging
10 - 1600 - 090	Skin Temperature Sensors 400 Series	90 / Case
11 - CA - 400 - 12	Reusable Temperature Adapter Cable (12 ft.)	Each



NEOTEMP

Disposable Neonatal Skin Temperature Sensors

		Part #	Packaging	Compatible With	Substitute For
GE/Ohmeda/Ohio		10 - 1930 - 020	20 / Box	Ohmeda Care Plus, I.C., CG Incubators & Infant Warming Systems 3000, 3300, 3500, 5000	OEM #LA003, LA005, #6600-0196-700 #6600-0208-700
		10 - 1930 - 060	60 / Case		
		10 - 1950 - 020	20 / Box	GE GIRAFFE OmniBed, GIRAFFE Incubators, GIRAFFE Warmers & PANDA Warmers	OEM #6600-0873-700 #6600-0874-700
		10 - 1950 - 060	60 / Case		
Dräger Medical/Air-Shields		10 - 1910 - 020	20 / Box	ISOLETTE Incubators (C-100, C-450, C-550) & PM-78 Radiant Warmers	OEM Premie Probe 1 #68-209-20/MU06933 #68-209-30/MU06935
		10 - 1910 - 060	60 / Case		
		10 - 1960 - 020	20 / Box	ISOLETTE C2000 & 8000 Incubators, VERSALET 7700 Care Center	OEM Premie Probe 5/5-L #83 -104 -15/MU12551
		10 - 1960 - 060	60 / Case		
		10 - 1965 - 020	20 / Box	CALEO Incubators & Babytherm Warmers 8004/8010	OEM #MX11000
	10 - 1965 - 060	60 / Case			
	10 - 1970 - 020	20 / Box	RESUSCITAIRE Radiant Warmers	OEM Premie Probe 3 #81-300-006, #81-300-008/MU11228 #81-300-009/MU11229	
	10 - 1970 - 060	60 / Case			
ATOM		10 - 2000D - 060	60 / Case	ATOM Medical Incubators & Infant Warmers	OEM NF-023, NF-024 #60882, #60883, #60084 w/ connector
		10 - 2000R	Packaged Individually		