

# **NOVA***lite*

## **LED Standard Single Use Stainless Steel Laryngoscope Blades and Handles**



- The ONLY LED Standard Single Use STAINLESS STEEL Laryngoscope Blades and Handles available
- Improved LED light output provides whiter, brighter light for enhanced visualization, ease of intubation
- LED is fixed to blade for improved patient safety
- LED is cool in operation – safer for neonates
- Complete range of blades: MacIntosh 0 - 5 and Miller 00 - 4
- Prevents cross contamination – improves reliability
- Eliminates costly reprocessing / cleaning / replacement

---

**NOVAMEDUSA**

**800-425-3535**  
**[www.NovaMed-USA.com](http://www.NovaMed-USA.com)**



# LED Standard Single Use Stainless Steel Laryngoscope Blades/Handles



- LED Illumination** Provides a whiter, brighter light for optimal illumination. Delivers 5X the light output of standard incandescent blades.
- 100% Stainless** The ONLY Standard 100% Stainless Steel Single Use LED Blades and Single Use Laryngoscope Handles available. 100% Stainless steel blade smoothly releases from handle, with minimal effort.
- Improved Safety** Addresses Joint Commission laryngoscope safety requirements. Single use blades / handles prevent cross-contamination / Improves reliability. Fixed light source - LED is fixed to blade for improved patient safety. LED is cool in operation – safer for neonates.
- Cost Effective** Single use eliminates reprocessing/cleaning costs. LED utilizes minimal energy consumption, saving on battery replacement. Blades are compatible with standard conventional handles.

MACINTOSH LED STANDARD SINGLE USE STAINLESS STEEL BLADES			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 7010D - 010	MAC 0, Neonate	55	9
18 - 7011D - 010	MAC 1, Infant	70	11
18 - 7012D - 010	MAC 2, Child	78	13
18 - 7013D - 010	MAC 3, Medium Adult	108	15
18 - 7013.5D - 010	MAC 3.5, Med Adult +	115	15
18 - 7014D - 010	MAC 4, Adult	133	15
18 - 7015D - 010	MAC 5, Large Adult	150	15

MILLER LED STANDARD SINGLE USE STAINLESS STEEL BLADES			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 70200D - 010	MIL 00, Premature	45	9.5
18 - 7020D - 010	MIL 0, Neonate	55	9.5
18 - 7021D - 010	MIL 1, Infant	80	9.5
18 - 7021.5D - 010	MIL 1.5, Infant +	90	9.5
18 - 7022D - 010	MIL 2, Child	130	13
18 - 7023D - 010	MIL 3, Medium Adult	170	13
18 - 7024D - 010	MIL 4, Adult	180	13

**LED Standard Single Use Laryngoscope Blades and Handles are packaged 10 pieces per box.**

Standard Laryngoscope Handles			
Order #	Description	Batteries	Packaging
18 - DMH - 010	Standard Single Use Medium/Adult Handle	2 C	10/Box
18 - DPH - 010	Standard Single Use Penlight/Pediatric Handle	2AAA	10/Box
14 - 600S - 01	Standard Reusable Medium/Adult Handle	2 C	Each
14 - 600S - 02	Standard Reusable Penlight/Pediatric Handle	2 AA	Each
14 - 600S - 03	Standard Reusable Stubby Handle	2 AA	Each