




LED Standard Laryngoscopes

The NEW Standard for Improved Illumination

- The ONLY LED SS blades available in Reusable & Single Use
- Whiter, brighter light for optimal task illumination
- 5X light output of standard incandescent blades
- Cool to the touch and low profile, ideal for neonatal



LED technology ensures whiter, brighter lighting with fewer shadows for optimal, focused task illumination. LED's offer significant clinical advantages with cost savings benefits.

- Reduced energy consumption prolongs battery life
 - Compatible with Standard Conventional handles
 - Upgrade to improved illumination at a cost savings
-

LED BRITE Reusable Standard SS Laryngoscope Blades



- Superior Illumination** Engineered with a whiter, brighter light.
5X light output for optimal, focused task illumination.
Ensures outstanding light transmission and improved visualization.
- Compatibility** Compatible with Standard Conventional Handles.
- Improved Safety** LED remains cool at all times.
Directional light is optimal for task performance.
Satin finish blades ensure ease of intubations to enhance patient safety.
- Cost Effective** No light filaments are utilized, product efficiency is optimized.
Reduced energy consumption minimizes battery replacement.
High quality stainless steel ensures long lasting performance.
- Comprehensive Selection** Available in both MacIntosh & Miller styles in a complete range of sizes.

MACINTOSH LED BRITE STANDARD BLADES			
Reusable Stainless Steel, Satin Finish			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 6010R	MAC 0, Neonate	55	9
18 - 6011R	MAC 1, Infant	70	11
18 - 6012R	MAC 2, Child	78	13
18 - 6013R	MAC 3, Medium Adult	108	15
18 - 6013.5R	MAC 3.5, Med Adult	115	15
18 - 6014R	MAC 4, Adult	133	15
18 - 6015R	MAC 5, Large Adult	150	15

MILLER LED BRITE STANDARD BLADES			
Reusable Stainless Steel, Satin Finish			
Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 60200R	MIL 00, Premature	45	9.5
18 - 6020R	MIL 0, Neonate	55	9.5
18 - 6021R	MIL 1, Infant	80	9.5
18 - 6021.5R	MIL 1.5, Infant	90	9.5
18 - 6022R	MIL 2, Child	130	13
18 - 6023R	MIL 3, Medium Adult	170	13
18 - 6024R	MIL 4, Adult	180	13

LED BRITE Single Use Standard SS Laryngoscope Blades



Whiter Brighter Light Engineered with a whiter, brighter light.
Improved directional, focused light with 5X light output.
Superior intubation conditions than with incandescent light sources.

Compatibility Compatible with Standard Conventional Handles.

Safe Prevents cross-infection.
Ideal for ER's or crash carts.
Individually wrapped, ready for use.
LED bulb remains cool at all times.
Directional light improves ease of intubations and patient safety.

Cost Effective Eliminates reprocessing/cleaning costs.
Reduced energy consumption minimizes battery replacement.

Comprehensive Selection Available in both MacIntosh & Miller styles in a complete range of sizes.

MACINTOSH LED BRITE STANDARD BLADES

Single Use Stainless Steel, Satin Finish

Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 7010D-010	MAC 0, Neonate	55	9
18 - 7011D-010	MAC 1, Infant	70	11
18 - 7012D-010	MAC 2, Child	78	13
18 - 7013D-010	MAC 3, Medium Adult	108	15
18 - 7013.5D-010	MAC 3.5, Med Adult	115	15
18 - 7014D-010	MAC 4, Adult	133	15
18 - 7015D-010	MAC 5, Large Adult	150	15

MILLER LED BRITE STANDARD BLADES

Single Use Stainless Steel, Satin Finish

Order #	Description/ Size	Length (mm)	Distal Width (mm)
18 - 7020D-010	MIL 00, Premature	45	9.5
18 - 7020D-010	MIL 0, Neonate	55	9.5
18 - 7021D-010	MIL 1, Infant	80	9.5
18 - 7021.5D-010	MIL 1.5, Infant	90	9.5
18 - 7022D-010	MIL 2, Child	130	13
18 - 7023D-010	MIL 3, Medium Adult	170	13
18 - 7024D-010	MIL 4, Adult	180	13

NOVALITE Standard Laryngoscopes Handles

Shown from left to right:
 Medium / Adult Handle
 Penlight / Pediatric Handle
 Stubby Handle
 Large Handle

Disposable Medium Handle shown below right.



- Reusable polished finish with high quality stainless steel construction – in four sizes to choose from.
- Compatible with other standard / conventional blades.
- Also available in a disposable stainless steel medium adult handle compatible with standard blades.



Order #	Description	Batteries
14 - 600S - 01	Reusable MEDIUM Adult Standard Conventional Laryngoscope Handle	2 C
14 - 600S - 02	Reusable PENLIGHT/ PEDIATRIC Standard Conventional Handle	2 AA
14 - 600S - 03	Reusable STUBBY Standard Conventional Laryngoscope Handle	2 AA
14 - 600S - 06	Reusable LARGE Standard Conventional Laryngoscope Handle	2 D
18-DMH-010	Single Use MEDIUM Adult Standard Laryngoscope Handle	2C
18-DPH-010	Single Use PENLIGHT Pediatric Standard Laryngoscope Handle	2AAA

Light Emitting Diode (LED) Replacement Lamps for Reusable LED Standard Laryngoscope Blades:

Order #	Description	Packaging
14 - 700LED - 06	Standard Reusable LED Replacement Lamp	6 per case