

NOVA*lite*



LED Standard Single Use Laryngoscope Blades and Handles



- The ONLY LED Standard Single Use STAINLESS STEEL Laryngoscope Blades available
- Improved LED light output provides whiter, brighter light for enhanced visualization, ease of intubation
- LED is fixed to blade for improved patient safety
- LED is cool in operation – safer for neonates
- Complete range of blades: MacIntosh 0 - 5 and Miller 00 - 4
- Prevents cross contamination – improves reliability
- Eliminates costly reprocessing / cleaning / replacement

NOVAMEDUSA

914-789-2100 / 800-425-3535
www.NovaMed-USA.com



LED Illumination

Provides a whiter, brighter light for optimal illumination. Delivers 5X the light output of standard incandescent blades.

100% Stainless

The ONLY Standard 100% Stainless Steel Single Use LED Blades available. Smoothly releases from handle, with minimal effort.

Improved Safety

Single use blades / handles prevent cross-contamination, improve reliability. Fixed light source - LED is fixed to blade for improved patient safety. LED is cool in operation – safer for neonates.

Cost Effective

Single use eliminates reprocessing/cleaning costs. LED utilizes minimal energy consumption, saving on battery replacement. Compatible with standard conventional handles.

MACINTOSH LED STANDARD SINGLE USE STAINLESS STEEL BLADES		
Order #	Description/ Size	Length (mm)
18 - 7010D - 010	MAC 0, Neonate	55
18 - 7011D - 010	MAC 1, Infant	70
18 - 7012D - 010	MAC 2, Child	78
18 - 7013D - 010	MAC 3, Medium Adult	108
18 - 7013.5D - 010	MAC 3.5, Med Adult +	115
18 - 7014D - 010	MAC 4, Adult	133
18 - 7015D - 010	MAC 5, Large Adult	150

MILLER LED STANDARD SINGLE USE STAINLESS STEEL BLADES		
Order #	Description/ Size	Length (mm)
18 - 70200D - 010	MIL 00, Premature	45
18 - 7020D - 010	MIL 0, Neonate	55
18 - 7021D - 010	MIL 1, Infant	80
18 - 7021.5D - 010	MIL 1.5, Infant +	90
18 - 7022D - 010	MIL 2, Child	130
18 - 7023D - 010	MIL 3, Medium Adult	170
18 - 7024D - 010	MIL 4, Adult	180

LED Standard Single Use Laryngoscope Blades and Handles are packaged 10 pieces per box.

Standard Laryngoscope Handles			
Order #	Description	Batteries	Packaging
18 - DMH - 010	Single Use Standard Medium/ Adult Handle	2 C	10/Box
18 - DPH - 010	Single Use Standard Penlight/ Pediatric Handle	2 AA	10/Box
14 - 600S - 01	Reusable Standard Medium/ Adult Handle	2 C	Each
14 - 600S - 02	Reusable Standard Penlight/ Pediatric Handle	2 AA	Each
14 - 600S - 03	Reusable Standard Stubby Handle	2 AA	Each