

NOVALITE

Fiber Optic Laryngoscopes

XENON Technology - Green System



All NOVALITE Fiber Optic Blades are
Certified Safe for the MRI Suite paired
with a NOVALITE MRI Conditional Handle.

NOVAMEDUSA

800-425-3535
www.NovaMed-USA.com

NOVALITE Fiber Optic Laryngoscopes

Available in Reusable and Single Use Fiber Optics with Ultra-Brite XENON Technology



► Illumination

Engineered with XENON light technology to facilitate superior illumination. Designed with an enlarged 4mm fiber optic bundle for improved light output.

► Compatibility

Utilizes the universal green system to ensure interchangeability; compatible with all green system handles and blades from competitive brands.

► Patient Safety

Satin finish blades enhance the ease of intubations and improve patient safety. Ensures enhanced patient comfort.

► Flexibility

Designed with an adjustable and entirely removable Fiber Optic bundle for enhanced cleaning benefits for improved safety and enhanced hospital risk management.

► Extensive Selection

Available in a comprehensive range of blade styles and sizes including MacIntosh 3.5, Miller 1.5, and Wis-Hipple 1.5.

Fiber Optic Handles



Xenon Light Technology Replacement Bulbs



Size	Batteries	Model #
Medium/Adult	2C	14 - 600FO - 01
Penlight/Pediatric	2AA	14 - 600FO - 02
Stubby	2AA	14 - 600FO - 03

Replacement Xenon Bulbs	Model #
Ultra Brite Xenon Bulbs 2.7V (6 pack)	14 - 700FO - 02 - 06
MRI Conditional 3.5V Xenon Bulbs (6 pack)	14 - 700FO - 03 - 06



MRI Conditional Handle / Battery	Model #
MRI Conditional Fiber Optic Medium Handle	15-MRI-HMED-1
Lithium Battery (MRI Conditional)	15-MRI-L1



NOVALITE Fiber Optic Reusable Laryngoscope Blades

Stainless Steel, Satin Finish – Green System

High quality stainless steel ensures long lasting performance.

Featuring an enlarged fiber optic bundle for optimal illumination.

MACINTOSH FIBER OPTIC REUSABLE

Model #	Size / Description
14 - 100FO - MAC - 0	0, Neonate
14 - 100FO - MAC - 1	1, Infant
14 - 100FO - MAC - 2	2, Child
14 - 100FO - MAC - 3	3, Medium Adult
14 - 100FO - MAC - 3.5	3.5, Medium Adult +
14 - 100FO - MAC - 4	4, Adult
14 - 100FO - MAC - 5	5, Large Adult

MILLER FIBER OPTIC REUSABLE

Model #	Size / Description
14 - 300FO - MIL - 00	00, Premature
14 - 300FO - MIL - 0	0, Neonate
14 - 300FO - MIL - 1	1, Infant
14 - 300FO - MIL - 1.5	1.5, Infant +
14 - 300FO - MIL - 2	2, Child
14 - 300FO - MIL - 3	3, Medium Adult
14 - 300FO - MIL - 4	4, Adult

FLIP TIP MACINTOSH FIBER OPTIC REUSABLE

Designed to facilitate
difficult intubations with
flexible tip adjustable
through 70°

Model #	Size / Description
14 - 200FOFT - MAC - 2	2, Child
14 - 200FOFT - MAC - 3	3, Medium Adult
14 - 200FOFT - MAC - 4	4, Adult

SPECIALTY BLADES FIBER OPTIC REUSABLE

	Model #	Size
Robert Shaw	14 - 400FO - RS - 0	0, Neonate
	14 - 400FO - RS - 1	1, Infant
Wis - Hipple	14 - 800FO - WH - 1.5	1.5, Infant
Cranwall	14 - 400FO - CW - 2	2, Child
	14 - 400FO - CW - 3	3, Medium Adult
	14 - 400FO - CW - 4	4, Adult

NOVALITE Disposable / Single Use Fiber Optic Options

- For safer intubation preventing against patient cross infections.
- Affordably priced to eliminate cleaning requirements.



Fiber Optic Low Profile Hybrid Laryngoscope Blades

- Single use Stainless Steel blade with polished Acrylic bundle to maximize illumination.
- Ideal for crash carts / intubation boxes, NICU, PICU and ER.
- Compatible with Green System Fiber Optic Handles.
- Packaged individually.

MACINTOSH FIBER OPTIC HYBRID	
Model #	Size / Description
16 - 100FOMD - MAC - 0	0, Neonate
16 - 100FOMD - MAC - 1	1, Infant
16 - 100FOMD - MAC - 2	2, Child
16 - 100FOMD - MAC - 3	3, Medium Adult
16 - 100FOMD - MAC - 4	4, Adult
16 - 100FOMD - MAC - 5	5, Large Adult

MILLER FIBER OPTIC HYBRID	
Model #	Size / Description
16 - 300FOMD - MIL - 00	00, Premature
16 - 300FOMD - MIL - 0	0, Neonate
16 - 300FOMD - MIL - 1	1, Infant
16 - 300FOMD - MIL - 2	2, Child
16 - 300FOMD - MIL - 3	3, Medium Adult
16 - 300FOMD - MIL - 4	4, Adult



Fiber Optic Disposable Laryngoscope Blades

- High quality plexi construction, with desired firmness for safer use.
- Polished stem Fiber Optic bundle for improved illumination.
- Compatible with Green System Fiber Optic Handles.
- Packaged 10 / pre-pack.

MACINTOSH FIBER OPTIC DISPOSABLE	
Model #	Size / Description
14 - 100FOD - MAC - 2	2, Child
14 - 100FOD - MAC - 3	3, Medium Adult
14 - 100FOD - MAC - 4	4, Adult

MILLER FIBER OPTIC DISPOSABLE	
Model #	Size / Description
14 - 300FOD - MIL - 0	0, Neonate
14 - 300FOD - MIL - 1	1, Infant
14 - 300FOD - MIL - 2	2, Child
14 - 300FOD - MIL - 3	3, Medium Adult