



## Fiber Optic Single Use Hybrid Laryngoscope Blade Stainless Steel Green System with Acrylic FO Bundle



- ✓ For safer intubation eliminating cross contamination.
- ✓ Affordably priced to eliminate cleaning requirements.
- ✓ Single use stainless steel blade with polished Acrylic bundle engineered to maximize illumination.
- ✓ Ideal for crash carts / intubation boxes, NICU, PICU, L&D and ER.
- ✓ Compatible with Green System Fiber Optic Handles.

<b>MACINTOSH FIBER OPTIC SINGLE USE</b>	
<b>Part # (10 pcs / Package)</b>	<b>Size / Description</b>
16 - 100FOMD - MAC - 0 - 010	0, Neonate
16 - 100FOMD - MAC - 1 - 010	1, Infant
16 - 100FOMD - MAC - 2 - 010	2, Child
16 - 100FOMD - MAC - 3 - 010	3, Medium Adult
16 - 100FOMD - MAC - 4 - 010	4, Adult
16 - 100FOMD - MAC - 5 - 010	5, Large Adult
<b>OXFORD FIBER OPTIC SINGLE USE</b>	
10 - 400FOMD - OX - 010	1, Infant

<b>MILLER FIBER OPTIC SINGLE USE</b>	
<b>Part # (10 pcs / Package)</b>	<b>Size / Description</b>
16 - 300FOMD - MIL - 00 - 010	00, Premie
16 - 300FOMD - MIL - 0 - 010	0, Neonate
16 - 300FOMD - MIL - 1 - 010	1, Infant
16 - 300FOMD - MIL - 2 - 010	2, Child
16 - 300FOMD - MIL - 3 - 010	3, Medium Adult
16 - 300FOMD - MIL - 4 - 010	4, Adult