



NOVAMEDUSA

Precise Temperature Monitoring for **AIRBORNE**
Life Support Systems

NOVATEMP Skin Temperature Sensors 400 Series

- Accuracy** NOVAMED thermistor sensor 400 Series guarantees continuous temperature reading sensitivity with accuracy $\pm 0.2^{\circ}\text{C}$ for the full scale range ($25^{\circ}\text{C} - 50^{\circ}\text{C}$).
- Safety** Latex-Free hypoallergenic one inch diameter adhesive foam disk is dermatologically and allergy-tested for even the most sensitive skin.
- Versatility** Ideal for baby-sensitive temperature monitoring.
- Quick Response** Unique thermistor design ensures rapid response to reflect any temperature change on transport patients.
- Stress-Free** NOVATEMP sensors are engineered with a longer lead wire and a custom fitted molded connector to further ensure response time reliability.

NOVATEMP SKIN TEMPERATURE SENSORS / 400 SERIES		
Order Ref #	Product	Packaging
10-1600-030	Skin Temperature Sensors/ 400 Series	30/ Box
10-1600-090	Skin Temperature Sensors/ 400 Series	90/ Case

NOVATEMP Temperature Adapter Cable for AirBorne Transport Equipment



Reusable Temperature Adapter Cable with 1/4" Right Angle Phono Plug Connector for AirBorne Transport Equipment and Molex Connector for 400 Series Temperature Sensors.

Available in 12' Length

Order Ref # 11-CA-400-12