



NOVAMEDUSA

The Preferred Source in Temperature Monitoring

NOVATEMP General Purpose Temperature Probes

User Friendly

Flexible, Latex-Free tubing enhances ease of insertion, promoting maximum patient comfort.

Complete Size Variety

Available in both 9FR and 12FR for a variety of clinical applications. NOVAMED recommends use for nasopharyngeal, esophageal, axillary or rectal placement.

Thermistor Accuracy

NOVAMED thermistor sensor - 400 Series – guarantees temperature reading sensitivity with accuracy +/- 0.2° C, full scale (25° - 50° C).

Stress-Free Reliability

Double strain relief, extended lead wire and custom fitted molded connector ensure operating room reliability.

NOVATEMP General Purpose Temperature Probes

Packaged 40 / box; 160 / case. Model numbers below reflect 160 per case.



400 Series	
Order Ref #	Size
10 - 1610 - 160	9FR
10 - 1620 - 160	12FR

NOVATEMP Temperature Adapter Cables

NovaMed offers a wide range of temperature adapter cables to match your needs according to the physiological monitors utilized.

Temperature Adapter Cable w/ Molex Connector, Terminating in:	Length	Order Ref. #
¼" Right Angle Phono Plug, 400 Series Thermistor	12'	11-CA-400-12
2 Pin Din Connector, use w/ Philips (Hewlett-Packard) monitors	12'	11-CA-401-12
¼" Straight Phono Plug, 400 Series Thermistor	12'	11-CA-402-12
7 Pin Din Connector, use w/ Dräger Medical/Siemens SC Series monitors	12'	11-CA-407-12